

Cyan fluorescent protein TagCFP

- Bright cyan fluorescence
- Monomer with successful performance in fusions
- High pH stability and photostability
- Proven suitability to generate stably transfected cell lines
- Recommended for protein labeling

Description

TagCFP is a cyan monomeric protein generated on the basis of GFP-like protein from jellyfish *Aequorea macrodactyla* (Xia et al., 2002). It possesses bright fluorescence with excitation/emission maxima at 458 and 480 nm, respectively. TagCFP is significantly brighter than commonly used ECFP. TagCFP is mainly intended for protein labeling in protein localization and interaction studies. It can also be used for cell and organelle labeling and for tracking the promoter activity.

Main properties of TagCFP

Characteristic	
Molecular weight	27 kDa
Polypeptide length	239 aa
Fluorescence color	cyan
Excitation max	458 nm
Emission max	480 nm
Quantum yield	0.57
Extinction coefficient	37 000 M ⁻¹ cm ⁻¹
Brightness*	21.1
Brightness % of EGFP	64
pKa	4.7
Structure	monomer
Aggregation	no
Maturation rate at 37°C	fast
Photostability	high

*Brightness is a product of extinction coefficient and quantum yield, divided by 1000.


Performance and use

TagCFP can be easily expressed and detected in a wide range of organisms. Mammalian cells transiently transfected with TagCFP expression vectors give bright fluorescent signals in 10-12 hrs after transfection. No cell toxic effects and visible protein aggregation are observed.

TagCFP performance in fusions has been demonstrated in human cytoplasmic beta-actin and alpha-tubulin models.


TagCFP suitability to generate stably transfected cells has been proven by Marinpharm company (www.marinpharm.com). Cell line expressing TagCFP fusion with mitochondrial targeting sequence (MTS) is commercially available.

TagCFP can be used in multicolor labeling applications with green, yellow, red, and far-red fluorescent dyes.



TagCFP normalized excitation (thin line) and emission (thick line) spectra.

Complete TagCFP spectra in Excel format can be downloaded from the Evrogen Web site at www.evrogen.com/products/TagCFP/TagCFP.shtml



Multicolor labeling of mammalian cells.

A — mitochondria-targeted TagCFP (cyan), Dendra2-vimentin fusion (green), and TagRFP-actin fusion (red); B — mitochondria-targeted TagCFP (cyan), peroxisome-targeted PhiYFP-m (yellow), and membrane-targeted KillerRed (red); C — TagCFP-actin fusion (blue), TagYFP-tubulin fusion (green), TTagFP635-H2B fusion (red), and Golgi-targeted TagRFP (orange).

Available variants and fusions

TagCFP codon usage is optimized for high expression in mammalian cells (Haas *et al.*, 1996), but it can be successfully expressed in many other heterologous systems.

TagCFP-mito fusion

A mitochondrial targeting sequence (MTS) is linked to the TagCFP N-terminus. MTS was derived from the subunit VIII of human cytochrome C oxidase (Rizzuto *et al.*, 1989; Rizzuto *et al.*, 1995). When expressed in mammalian cells, this variant provides cyan fluorescent labeling of mitochondria.

TagCFP-actin fusion

Human beta-actin is linked to the TagCFP C-terminus. When expressed in mammalian cells, this fusion provides cyan fluorescent labeling of actin filaments.

TagCFP-tubulin fusion

Human alpha-tubulin is linked to the TagCFP C-terminus. When expressed in mammalian cells, this fusion provides cyan fluorescent labeling of tubulin filaments.

Recommended filter sets and antibodies

TagCFP can be detected using fluorescence filter sets for ECFP and the similar. Recommended Omega Optical filter sets are XF114-2 and XF130-2.

TagCFP can be recognized using Anti-Tag(CGY)FP antibody (Cat.# AB121-AB122) available from Evrogen.




TagCFP licensing opportunities

Evrogen technology embodied in TagCFP is available for expanded and commercial use with an adaptable licensing program. Benefits from flexible and market-driven license options are offered for upgrade and novel development of products and applications. For licensing information, please contact Evrogen at license@evrogen.com.

TagCFP comprises the following amino acid substitutions compared with wild-type *A. macrodactyla* GFP (AY013824): K3G, T9A, F64L, S65A, Y66W, F99H, I123V, M128E, D129G, N144S, F145A, N146I, H148D, K162E, V163A, T203C, T205S, T214A, F220L, F223S, C227Y, G228C, K238R. It has 77% amino acid sequence identity with wild-type GFP from *A. victoria*.


References

- Haas *et al.* (1996) *Curr. Biol.* 6: 315–324.
- Rizzuto *et al.* (1989) *J. Biol. Chem.* 264: 10595–10600.
- Rizzuto *et al.* (1995) *Curr. Biol.* 5: 635–642.
- Xia *et al.* (2002) *Mar. Biotechnol.* 4(2): 155-162.



Expression of TagCFP fusions in transiently transfected mammalian cells (HeLa).

Confocal microscopy: A — mitochondria-targeted TagCFP; B — TagCFP fusion with alpha-tubulin; C — TagCFP fusion with cytoplasmic beta-actin.



Expression of mitochondria-targeted TagCFP in stably transfected U-205 human osteosarcoma cells.

Photograph of stably transfected cell line was provided by Dr. Christian Petzelt (Marinpharm).

TagCFP-related products

Product	Cat.#	Description	Size
TagCFP expression/source vectors			
pTagCFP-C	FP111	Mammalian expression vector encoding humanized TagCFP and allowing TagCFP expression and generation of fusions to the TagCFP C-terminus	20 µg
pTagCFP-N	FP112	Mammalian expression vector encoding humanized TagCFP and allowing TagCFP expression and generation of fusions to the TagCFP N-terminus	20 µg
pTagCFP-actin	FP114	Mammalian expression vector encoding humanized TagCFP fused with human cytoplasmic beta-actin	20 µg
pTagCFP-tubulin	FP115	Mammalian expression vector encoding humanized TagCFP fused with human alpha-tubulin	20 µg
pTagCFP-mito	FP117	Mammalian expression vector encoding humanized TagCFP fused with mitochondria localization signal	20 µg
pTagCFP-keratin	FP119	Mammalian expression vector encoding humanized TagCFP fused with human cytokeratin-18	20 µg
Antibodies against TagCFP			
Anti-Tag(CGY)FP antibody	AB121	Rabbit polyclonal antibody against TagCFP, TagGFP,	100 µg
	AB122	TagGFP2, TagYFP, PS-CFP2, Case12, HyPer, and EGFP	200 µg

Please contact your local distributor for exact prices and delivery information.

Notice to Purchaser:

TagCFP-related products: These products are intended for research use only and covered by Evrogen Patents and/or Patent applications pending. By use of these products, you accept the terms and conditions of the applicable Limited Use Label License (available at www.evrogen.com//Evrogen-FP-license.shtml).

CMV Promoter: The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839, and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242.