


Bacterial expression vector pTurboYFP-B

Product	Cat.#	Size
pTurboYFP-B	FP613	20 µg

Please contact your local distributor for exact prices and delivery information.



Vector type	bacterial expression vector
Reporter	TurboYFP
Reporter codon usage	mammalian
Promoter for TurboYFP	T5 promoter/lac operator
Host cells	prokaryotic
Selection	ampicillin
Replication	ColE1 ori

For vector sequence, please visit our Web site at www.evrogen.com/support/vector-info.shtml

Use

- TurboYFP expression in bacterial cells using T5 promoter/lac operator
- Source of the TurboYFP coding sequence

References

Haas, J., et al. (1996) Codon usage limitation in the expression of HIV-1 envelope glycoprotein. *Curr. Biol.* 6:315–324.

Vector description

pTurboYFP-B is a bacterial expression vector encoding yellow fluorescent protein TurboYFP. Reporter codon usage is optimized for high expression in mammalian cells (humanized) (Haas et al., 1996).

The vector is primarily intended as a source of TurboYFP coding sequence. Flanking restriction sites are convenient for TurboYFP gene excision and its further insertion into other expression vectors of choice. Alternatively, TurboYFP coding sequence can be amplified by PCR.

Note: The plasmid DNA was isolated from dam⁺-methylated *E.coli*. Therefore some restriction sites are blocked by methylation. If you wish to digest the vector using such sites you will need to transform the vector into a dam⁻ host and make fresh DNA.

The vector can be also used for TurboYFP expression in prokaryotes under the control of T5 promoter/lac operator. The vector backbone contains ColE1 origin of replication and ampicillin resistance gene for propagation and selection in *E. coli*.

Location of features:

T5 promoter/lac operator element: 7-87

T5 transcription start: 61

TurboYFP coding sequence: 133-864

Lambda t0 transcriptional termination region: 886-980

rrnB T1 transcriptional termination region: 1742-1840

ColE1 origin of replication: 2316

beta-lactamase coding sequence: 3934-3074

Notice to Purchaser:

Evrogen Fluorescent Protein Products (the Products) are intended for research use only and covered by Evrogen Patents and/or Patent applications pending. By use of these products, you accept the terms and conditions of the applicable Limited Use Label License (enclosed).

MATERIAL SAFETY DATA SHEET INFORMATION

To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.