


pTagCFP-tubulin vector

The vector sequence has been compiled using the information from sequence databases, published literature, and other sources, together with partial sequences obtained by Evrogen. This vector has not been completely sequenced.



For vector sequence, please visit our Web site at <http://www.evrogen.com/support/vector-info.shtml>

Location of features

PCMV IE: 1-589
 Enhancer region: 59-465
 TATA box: 554-560
 Transcription start point: 583
 TagCFP
 Kozak consensus translation initiation site: 606-616
 Start codon (ATG): 613-615; Stop codon: 2698-2700
 Last amino acid in TagCFP: 1324-1326
 Tubulin: 1345-2700
 SV40 early mRNA polyadenylation signal
 Polyadenylation signals: 2861-2866 & 2890-2895
 mRNA 3' ends: 2899 & 2911
 f1 single-strand DNA origin: 2958-3413
 Eukaryotic promoter for expression of Kan^r gene
 -35 region: 3475-3480; -10 region: 3498-3503
 Transcription start point: 3510
 SV40 origin of replication: 3754-3889
 SV40 early promoter
 Enhancer (72-bp tandem repeats): 3587-3658 & 3659-3730
 21-bp repeats: 3734-3754, 3755-3775 & 3777-3797
 Early promoter element: 3810-3816
 Major transcription start points: 3806, 3844, 3850 & 3855
 Kanamycin/neomycin resistance gene
 Neomycin phosphotransferase coding sequences:
 Start codon (ATG): 3938-3940; Stop codon: 4730-4732
 G->A mutation to remove Pst I site: 4120
 C->A (Arg to Ser) mutation to remove BssH II site: 4466
 Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal
 Polyadenylation signals: 4968-4973 & 4981-4986
 pUC plasmid replication origin: 5317-5960

Product	Cat.#	Size
pTagCFP-tubulin vector	FP115	20 µg
The price does not include delivery. The price varies in different countries. Please contact your local distributor for exact prices and delivery information.		
Vector type	mammalian expression vector	
Reporter	TagCFP	
Reporter codon usage	mammalian	
Promoter for TagCFP	P _{CMV IE}	
Host cells	mammalian	
Selection	prokaryotic - kanamycin eukaryotic - neomycin (G418)	
Replication	prokaryotic - pUC ori eukaryotic - SV40 ori	
Use	cyan fluorescent labeling of α-tubulin filaments	

Vector description

pTagCFP-tubulin is a mammalian expression vector encoding TagCFP-tubulin fusion protein. The vector can be used for fluorescent labeling of α-tubulin in living cells.

TagCFP codon usage is optimized for high expression in mammalian cells, i.e. humanized [Haas et al. 1996]. Human α-tubulin is fused to the TagCFP C-terminus. To increase mRNA translation efficiency, Kozak consensus translation initiation site is generated upstream of TagCFP-tubulin coding sequence [Kozak 1987].

pTagCFP-tubulin can be used as a source of TagCFP-tubulin hybrid sequence. The vector backbone contains unique restriction sites that permit its excision and further insertion into expression vector of choice.

Note: The plasmid DNA was isolated from dam⁺-methylated *E.coli*. Therefore some restriction sites are blocked by methylation. If you wish to digest the vector using such sites you will need to transform the vector into a dam⁻ host and make fresh DNA.

The vector backbone also contains immediate early promoter of cytomegalovirus (P_{CMV IE}) for protein expression, SV40 origin for replication in mammalian cells expressing SV40 T-antigen, pUC origin of replication for propagation in *E. coli*, and f1 origin for single-stranded DNA production. SV40 polyadenylation signals (SV40 poly A) direct proper processing of the 3' end of the reporter mRNA.

SV40 early promoter (P_{SV40}) provides neomycin resistance gene (Neo^r) expression to select stably transfected eukaryotic cells using G418. Bacterial promoter (P) provides kanamycin resistance gene expression (Kan^r) in *E. coli*. Kan^r/Neo^r gene is linked with herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signals.

Expression in mammalian cells

pTagCFP-tubulin can be transfected into mammalian cells by any known transfection method. CMV promoter provides strong, constitutive expression of the TagCFP-tubulin fusion in eukaryotic cells. If required, stable transformants can be selected using G418 [Gorman 1985].

Propagation in *E. coli*

Suitable host strains for propagation in *E. coli* include DH5alpha, HB101, XL1-Blue, and other general purpose strains. Plasmid incompatibility group is pMB1/ColE1. The vector confers resistance to kanamycin (30 µg/ml) to *E. coli* hosts. Copy number in *E. coli* is about 500.

References

- Gorman (1985). "High efficiency gene transfer into mammalian cells." In: *DNA cloning: A Practical Approach, Vol. II*. Ed. by Glover. (IRL Press, Oxford, U.K.) Pp. 143-190.
- Haas et al. (1996) "Codon usage limitation in the expression of HIV-1 envelope glycoprotein." *Curr Biol*, 6 (3): 315-324 / pmid: 8805248
- Kozak (1987) "An analysis of 5'-noncoding sequences from 699 vertebrate messenger RNAs." *Nucleic Acids Res*, 15 (20): 8125-8148 / pmid: 3313277

Notice to Purchaser:

Evrogen Fluorescent Protein Products (the Products) are intended for research use only. The Products are covered by U.S. Pat. 7,417,131 and other Evrogen Patents and/or Patent applications pending. By use of these Products, you accept the terms and conditions of the applicable Limited Use Label License.

The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839, and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242.

MATERIAL SAFETY DATA SHEET INFORMATION: To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.